

Your Power for Health



Your Power for Health



For further information, please visit our website [www.gbo.com/preanalytics](http://www.gbo.com/preanalytics) or contact us:

**Austria (Main Office)**  
Greiner Bio-One GmbH  
Phone (+43) 75 83 67 91-0  
Fax (+43) 75 83 63 18  
E-Mail [office@at.gbo.com](mailto:office@at.gbo.com)

**Brazil**  
Greiner Bio-One Brasil  
Phone (+55) 19 34 68 96 00  
Fax (+55) 19 34 68 96 21  
E-Mail [office@br.gbo.com](mailto:office@br.gbo.com)

**Egypt**  
Greiner Bio-One Middle East  
Phone (+20) 26 21 87 06  
Fax (+20) 26 21 87 08  
E-Mail [hisham.ouda@gbo.com](mailto:hisham.ouda@gbo.com)

**France**  
Greiner Bio-One SAS  
Phone (+33) 1 69 86 25 50  
Fax (+33) 1 69 86 25 36  
E-Mail [office@fr.gbo.com](mailto:office@fr.gbo.com)

**Germany**  
Greiner Bio-One GmbH  
Phone (+49) 201 861 86-0  
Fax (+49) 201 861 86-12  
E-Mail [office@de.gbo.com](mailto:office@de.gbo.com)

**Hungary**  
Greiner Bio-One Hungary Kft.  
Phone (+36) 96 21 30 88  
Fax (+36) 96 21 31 98  
E-Mail [office@hu.gbo.com](mailto:office@hu.gbo.com)

**Netherlands**  
Greiner Bio-One B.V.  
Phone (+31) 1 72 42 09 00  
Fax (+31) 1 72 44 38 01  
E-Mail [info@nl.gbo.com](mailto:info@nl.gbo.com)

**Spain**  
VACUETTE Espana S. A.  
Phone (+34) 91 652 77 07  
Fax (+34) 91 652 33 35  
E-Mail [info@vacuette.es](mailto:info@vacuette.es)

**Switzerland**  
Greiner Bio-One VACUETTE Schweiz GmbH  
Phone (+41) 7 12 28 55 22  
Fax (+41) 7 12 28 55 21  
E-Mail [office@ch.gbo.com](mailto:office@ch.gbo.com)

**Thailand**  
Greiner Bio-One Thailand Ltd  
Phone (+66) 3 84 65 63-30  
Fax (+66) 3 84 65 63-6  
E-Mail [office@th.gbo.com](mailto:office@th.gbo.com)

**UK**  
Greiner Bio-One Ltd.  
Phone (+44) 14 53 82 52 55  
Fax (+44) 14 53 82 62 66  
E-Mail [info@uk.gbo.com](mailto:info@uk.gbo.com)

**USA**  
Greiner Bio-One North America Inc.  
Phone (+1) 70 42 61 78 00  
Fax (+1) 70 42 61 78 99  
E-Mail [info@us.gbo.com](mailto:info@us.gbo.com)

980022 Rev. 02.03.2009 e



# VACUETTE® PREMIUM Tubes

Our Innovations  
for Your Safety

# VACUETTE® PREMIUM Tube

## PREMIUM Safety and PREMIUM Comfort

### The PREMIUM Tube



Since blood is potentially infectious sample material, our aim at Greiner Bio-One is to make daily routine procedures as simple and safe as possible for medical staff.

This was the stimulus for our product experts to develop the new, innovative **VACUETTE® PREMIUM Tube** with Safety Twist Cap. This may sound simple, but from the production technology point of view, this is a significant advance for user safety, requiring a high level of innovative strength and competence.

Maximum Handling Comfort	
<b>Easy manual opening</b>	The unique screw thread ensures effortless opening of the tube.
<b>Easy automatic opening</b>	The new twist cap is perfect for automated tube decapping.
<b>Better optical control of the blood sample</b>	The transparent plastic label allows better visual control of the sample at first glance.*
<b>One-handed opening</b>	If necessary, the tube can be easily opened with one hand.

Leading the market with technology is the prerequisite for developing and manufacturing these outstanding products. Greiner Bio-One is the only manufacturer in the world to offer an evacuated blood collection tube with Safety Twist Cap. This is once again the proof that Greiner Bio-One is Trendsetter for safety solutions in the field of blood collection.


With this innovation, handling comfort and safety are united in one product. It is the next logical step for even more safety when processing blood samples for analysis.

The **VACUETTE®** blood collection tube is already available in the PREMIUM version for most areas of application. For more detailed information, please refer to the **VACUETTE®** Product Catalogue.

Maximum Safety	
<b>No danger of blood splashes</b>	The gentle, controlled opening movement reduces the danger of blood splashes.
<b>Transport safety</b>	The Safety Twist Cap ensures safe transport, as the cap will not become loose, whatever the circumstances.
<b>Minimises aerosol effect</b>	No aerosol when opening due to the controlled, consistent opening movement.

\* Optional: the **VACUETTE®** PREMIUM tube is available in a choice of plastic label or standard paper label with writing space

## PREMIUM Safety



## PREMIUM Comfort

