

Your **Power** for Health



Your **Power** for Health



For further information, please visit our website www.gbo.com/preanalytics or contact us:

Austria (Main Office)

Greiner Bio-One GmbH
Phone (+43) 75 83 67 91-0
Fax (+43) 75 83 63 18
E-Mail office@at.gbo.com

France

Greiner Bio-One SAS
Phone (+33) 1 69 86 25 25
Fax (+33) 1 69 86 25 92
E-Mail office@fr.gbo.com

Netherlands

Greiner Bio-One B.V.
Phone (+31) 1 72 42 09 00
Fax (+31) 1 72 44 38 01
E-Mail info@nl.gbo.com

Thailand

Greiner Bio-One Thailand Ltd
Phone (+66) 3 84 65 63-30
Fax (+66) 3 84 65 63-6
E-Mail office@th.gbo.com

Brazil

Greiner Bio-One Brasil
Phone (+55) 19 34 68 96 00
Fax (+55) 19 34 68 96 21
E-Mail office@br.gbo.com

Germany

Greiner Bio-One GmbH
Phone (+49) 201 861 86-0
Fax (+49) 201 861 86-12
E-Mail office@de.gbo.com

Spain

VACUETTE Espana S. A.
Phone (+34) 91 652 77 07
Fax (+34) 91 652 33 35
E-Mail info@vacuette.es

UK

Greiner Bio-One Ltd.
Phone (+44) 14 53 82 52 55
Fax (+44) 14 53 82 62 66
E-Mail info@uk.gbo.com

Egypt

Greiner Bio-One Middle East
Phone (+20) 26 21 87 06
Fax (+20) 26 21 87 08
E-Mail hisham.ouda@gbo.com

Hungary

Greiner Bio-One Hungary Kft.
Phone (+36) 96 21 30 88
Fax (+36) 96 21 31 98
E-Mail office@hu.gbo.com

Switzerland

Greiner Bio-One VACUETTE
Schweiz GmbH
Phone (+41) 7 12 28 55 22
Fax (+41) 7 12 28 55 21
E-Mail office@ch.gbo.com

USA

Greiner Bio-One North America Inc.
Phone (+1) 70 42 61 78 00
Fax (+1) 70 42 61 78 99
E-Mail info@us.gbo.com

980184_Rev01_12.10 e



VACUETTE® VISIO PLUS Needle

**Bringing Vision to
Blood Collection**

VACUETTE® VISIO PLUS Needle

The **VACUETTE® VISIO PLUS Needle** with optical venipuncture control can make an important contribution to simplifying blood collection routine for phlebotomists.

The needle has a translucent plastic hub, which functions as a view window.

If venipuncture is successful, the view window is coloured red due to the blood flow.

The VISIO PLUS Needle guarantees a safe blood collection procedure for the patient, whilst at the same time ensuring best possible comfort. Optical venipuncture control via the view window makes it much easier to find the vein lumen.

Re-positioning in the vein which may cause patient discomfort, can be avoided with the VISIO PLUS venipuncture aid.

Perfected grinding technology allows for sharp needle tip and optimal tissue penetration.

To be able to differentiate from the VACUETTE® standard needles, the VISIO PLUS blood collection needles have transparent needle containers.

VISIO PLUS – the Greiner Bio-One innovation.

Item n°	Description	Size	Packaging Box / Carton
450040	VACUETTE® VISIO PLUS Needle green, sterile	0.8 x 38 mm / 21G x 1 ½"	100 pcs./2000 pcs.
450041	VACUETTE® VISIO PLUS Needle black, sterile	0.7 x 38 mm / 22G x 1 ½"	100 pcs./2000 pcs.
450042	VACUETTE® VISIO PLUS Needle green, sterile	0.8 x 25 mm / 21G x 1"	100 pcs./2000 pcs.
450043	VACUETTE® VISIO PLUS Needle black, sterile	0.7 x 25 mm / 22G x 1"	100 pcs./2000 pcs.
450235	VACUETTE® QUICKSHIELD Complete PLUS, sterile, single-packed <small>Safety Tube Holder & VACUETTE® VISIO PLUS Needle Combination Product</small>	0.8 x 25 mm / 21G x 1"	30 pcs./300 pcs.
450239	VACUETTE® QUICKSHIELD Complete PLUS, sterile, single-packed <small>Safety Tube Holder & VACUETTE® VISIO PLUS Needle Combination Product</small>	0.8 x 38 mm / 21G x 1 ½"	30 pcs./300 pcs.
450226	VACUETTE® QUICKSHIELD Complete PLUS, sterile, single-packed <small>Safety Tube Holder & VACUETTE® VISIO PLUS Needle Combination Product</small>	0.7 x 25 mm / 22G x 1"	30 pcs./300 pcs.
450228	VACUETTE® QUICKSHIELD Complete PLUS, sterile, single-packed <small>Safety Tube Holder & VACUETTE® VISIO PLUS Needle Combination Product</small>	0.7 x 38 mm / 22G x 1 ½"	30 pcs./300 pcs.

VACUETTE® QUICKSHIELD Complete PLUS

Simply Practical.

The VACUETTE® QUICKSHIELD Safety Tube Holder is especially suitable for routine blood collection. This product can provide the user with the simplest handling and reliable infection protection with no change to usual collection technique. The pre-assembled product, safety tube holder with the VACUETTE® VISIO PLUS needle, has the advantage of the following features:

- The holder comes with pre-attached needle, sterile and single-packed.
- Transparent view window for optical venipuncture control.
- It is not necessary to thread in the needle.
- A practical, ready-to-use solution that saves time.
- Handling is optimised.

Safe disposal due to protected needle.

The view window made out of transparent plastic allows for optical control of the blood flow.

Perfected grinding technology allows for sharp needle tip and optimal tissue penetration.

